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Installation Instructions for Julian USA Ceiling Unit

General:

The Julian USA Ceiling Unit has been equipped with reserve cylinders in order to conform with standards. Due to packaging complications, the cylinder holder and the pressure reducer must be mounted during on-site installation. In order to ensure that assembly occurs correctly, i.e. that hoses for the pressure reducer are laid correctly (hoses and cables must be correctly secured to prevent kinks at the ceiling connection), the following steps must be taken during installation.

Required Accessories:

- 1 cylinder guides, complete
- O₂ pressure reducer
- N₂O pressure reducer (optional), if applicable
- Julian USA Ceiling Unit

Installation Instructions:

8. Mount both cylinder guides to Julian as show in figure 1. Use screws and washers to secure them to the castor receptacle and to secure them to the tie-bar.
Note: One cylinder holder has a long screw on the upper strap, see figure 2. This shall be mounted below the gas inlet block.
9. Mount pressure reducer(s) O₂ / (N₂O) to the gas inlet block and plug the connecting plug of the high-pressure sensor into the jacks at the back of the device according to the code.
10. If reserve cylinder(s) are present, mount cylinders to the pre-mounted cylinder guide. Place the adhesive safety strips on the cylinder holder to prevent scratching.
11. Secure the outlet hose and cable of the O₂ pressure reducer to the cylinder holder below the gas inlet block using the hose clamp as shown in figure 2. To do this, loosen the (long) screw on the strap and mount it along with the hose clamp. The hoses must be guided in such a way as to ensure that both hose and cable are not too close to the ceiling receptacle (the area between the cylinders).
12. Place the outlet hose and cable of the N₂O pressure reducer as shown in illustration 3 (if applicable). To do this, remove the guiding sheet from the tie-bar as shown in figure 3. Place the hose and cable in the tie-bar and secure using the guiding sheet. Now secure the whole (hose and cable) to the back of the monitor using the hose clamp, witch is mounted on the back wall.

The hoses must be guided in such a way as to ensure that both hose and cable are not too close to the ceiling receptacle (the area between the cylinders).
13. If cylinders are present, mount pressure reducer to the cylinders.
14. Give these installation instructions and the Operator's Manual to the customer.

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without express authority. All rights reserved.
Dräger Medizintechnik GmbH, D-23542 Lübeck
Dräger Medizintechnik GmbH, D-23542 Lubeck